


Anti-Slfn12 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Iym2

Item Number: CM0019689



Introduction

High-quality rabbit polyclonal antibody against human SLFN12, validated for Western blot and ELISA. Ideal for studying estrogen-induced apoptosis, rRNA degradation, and cell differentiation pathways. SWISS: Q8IYM2.

[Learn More](#)

Parameter	Value
Product Name	SLFN12 Rabbit pAb
Remarks/Alias	SLFN3; SLFN12
Species	Human
Gene ID (Human)	55106
Gene ID	55106
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-390 of human SLFN12 (NP_001275938.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IYM2
Protein Weight	67kDa
Shipping	Ice bag