


Anti-Slirp Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Gzt3

Item Number: CM0018519



Introduction

Rabbit polyclonal antibody against human SLIRP, a mitochondrial RNA-binding protein and nuclear receptor corepressor. Validated for WB, IHC-P, ELISA with cross-reactivity to human and mouse. Ideal for studying SRA-mediated transcriptional regulation and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	SLIRP Rabbit pAb
Remarks/Alias	DC50; PD04872; C14orf156; SLIRP
Species	Human
GeneID (Human)	81892
GeneID	81892
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLIRP (NP_112487.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9GZT3
Protein Weight	12kDa
Shipping	Ice bag