

Anti-Slit1 Rabbit Polyclonal Antibody For Wb, Elisa - O75093

Item Number: CM0022520

Introduction

SLIT1 Rabbit Polyclonal Antibody (O75093) validated for WB and ELISA. Recognizes human SLIT1, a secreted guidance cue involved in neural development. Ideal for neuroscience research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	SLIT1 Rabbit pAb
Remarks/Alias	MEGF4; SLIL1; SLIT3; SLIT-1; SLIT1
Species	Human
GeneID (human)	6585
GeneID	6585
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1315-1534 of human SLIT1 (NP_003052.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75093
Protein Weight	168kDa
Shipping	Ice bag