


Anti-Slit2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O94813

Item Number: CM0010573



Introduction

Recombinant rabbit monoclonal antibody against human Slit2 (O94813). Validated for WB, IHC-P, and ELISA with cross-reactivity to mouse and rat. Targets a secreted guidance cue involved in axonal navigation and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Slit2 Rabbit mAb
Remarks/Alias	SLIL3, Slit-2, Slit2
Species	Human
GeneID (Human)	9353
GeneID	9353
Immunogen	Synthetic peptide corresponding to amino acids 1400-1500 of human Slit2 (O94813)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94813
Protein Weight	170 kDa
Shipping	Ice bag