

Anti-Slmap Rabbit Polyclonal Antibody For Wb, Elisa - Q14Bn4

Item Number: CM0022018



Introduction

Rabbit polyclonal antibody against human SLMAP protein. Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Useful for investigating Hippo signaling and cell growth control.

[Learn More](#)

Parameter	Value
Product Name	SLMAP Rabbit pAb
Remarks/Alias	SLAP; SLMAP
Species	Human
Gene ID (Human)	7871
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 540-790 of human SLMAP (NP_009090.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14BN4
Protein Weight	95kDa
Shipping	Ice bag