

## Anti-Slp2 Rabbit Monoclonal Antibody For Western Blot And Elisa - Q9Ujz1

Item Number: CM0003237

### Introduction

High-specificity anti-SLP2 rabbit monoclonal antibody validated for WB and ELISA. Targets human SLP2 with cross-reactivity to mouse and rat. Ideal for mitochondrial research and IVD assay development.

[Learn More](#)

| Parameter        | Value                                                                                                       |
|------------------|-------------------------------------------------------------------------------------------------------------|
| Product Name     | SLP2 Rabbit mAb                                                                                             |
| Remarks/Alias    | SLP-2; HSPC108; SLP2                                                                                        |
| Species          | Human                                                                                                       |
| GeneID (Human)   | 30968                                                                                                       |
| GeneID           | 30968                                                                                                       |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 10-130 of human SLP2 (Q9UJZ1) |
| Source           | Rabbit                                                                                                      |
| Category         | Monoclonal Antibodies                                                                                       |
| Application      | WB, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                           |
| SWISS            | Q9UJZ1                                                                                                      |
| Protein Weight   | 39kDa                                                                                                       |
| Shipping         | Ice bag                                                                                                     |