

Anti-Slpi Rabbit Polyclonal Antibody For Wb, Elisa - P03973

Item Number: CM0024384



Introduction

SLPI rabbit polyclonal antibody validated for WB, ELISA. Detects antileukoproteinase in human, mouse, rat; a 14 kDa protease inhibitor of innate immunity. For IVD research.

[Learn More](#)

Parameter	Value
Product Name	SLPI Rabbit pAb
Remarks/Alias	ALP; MPI; ALK1; BLPI; HUSI; WAP4; WFDC4; HUSI-I; SLPI
Species	Human
GeneID (Human)	6590
GeneID	6590
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03973
Protein Weight	14kDa
Shipping	Ice bag