

Anti-Sltm Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nwh9

Item Number: CM0023871

Introduction

Rabbit polyclonal antibody targeting human SLTM (SAFB-like transcription modulator), validated for WB and ELISA. Cross-reacts with mouse. Suitable for transcription regulation and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	SLTM Rabbit pAb
Remarks/Alias	Met; SLTM
Species	Human
GeneID (Human)	79811
GeneID	79811
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 900-1000 of human SLTM (NP_079031.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NWH9
Protein Weight	117kDa
Shipping	Ice bag