

Anti-Slug Polyclonal Antibody For Wb And Elisa - O43623

Item Number: CM0015127

Introduction

Slug Rabbit pAb is a polyclonal antibody raised against recombinant human Slug (aa1-150). Validated for WB and ELISA, reacts with human and mouse. Ideal for studying EMT, neural crest development, and cancer metastasis.

[Learn More](#)

Parameter	Value
Product Name	Slug Rabbit pAb
Remarks/Alias	SLUG; WS2D; SLUGH; SLUGH1; SNAIL2; Slug
Species	Human
GeneID (human)	6591
GeneID	6591
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human Slug (NP_003059.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43623
Protein Weight	30kDa
Shipping	Ice bag