

Anti-Slug Rabbit Monoclonal Antibody For Wb And Elisa - O43623

Item Number: CM0010862

Introduction

Rabbit monoclonal antibody targeting human Slug (SNAI2) for WB and ELISA. Cross-reacts with human and mouse. Ideal for EMT, neural crest development, and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	Slug Rabbit mAb
Remarks/Alias	SLUG; WS2D; SLUGH; SLUGH1; SNAI2
Species	Human
GenelD (human)	6591
GenelD	6591
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43623
Protein Weight	30kDa
Shipping	Ice bag