

Anti-Slug Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43623

Item Number: CM0018297

Introduction

Anti-Slug Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Detects Slug/SNAI2 in Human, Mouse, Rat. Ideal for EMT and cancer research. - O43623

[Learn More](#)

Parameter	Value
Product Name	Slug Rabbit pAb
Remarks/Alias	SLUG; WS2D; SLUGH; SLUGH1; SNAI2; Slug
Species	Human
Gene ID (Human)	6591
Gene ID	6591
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 21-120 of human Slug (NP_003059.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43623
Protein Weight	30kDa
Shipping	Ice bag