


## Anti-Smad1 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q15797

Item Number: CM0003193



### Introduction

Rabbit monoclonal antibody against Smad1, validated for Western blot, IHC-P, IP, and ELISA. Recognizes human and mouse Smad1; suitable for BMP/TGF- $\beta$  signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] Smad1 Rabbit mAb                                                                                    |
| Remarks/Alias    | BSP1, JV41, BSP-1, JV4-1, MADH1, MADR1, d1                                                                         |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 4086                                                                                                               |
| GeneID           | 4086                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 20-240 of human Smad1 (NP_005891.1). |
| Source           | Rabbit                                                                                                             |
| Category         | Monoclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IP, ELISA                                                                                               |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | Q15797                                                                                                             |
| Protein Weight   | 52kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |