


Anti-Smad1 Monoclonal Antibody (Ko Validated) For Wb, If/Icc, Elisa - Q15797

Item Number: CM0007686



Introduction

Rabbit monoclonal antibody validated by knockout (KO) against human Smad1, suitable for Western blot (WB), immunofluorescence/immunocytochemistry (IF/ICC), and ELISA. Detects Smad1 in human and mouse species. Smad1 is a key transcriptional modulator in BMP signaling pathways involved in embryonic development and cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Smad1 Rabbit mAb
Remarks/Alias	BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1; d1
Species	Human
GenelD (Human)	4086
GenelD	4086
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 160-260 of human Smad1 (NP_001003688.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15797
Protein Weight	52kDa
Shipping	Ice bag