


Anti-Smad1 Rabbit Monoclonal Antibody (Ko Validated) For Wb, Ihc-P, Ip, Elisa - Q15797

Item Number: CM0006807



Introduction

KO-validated rabbit monoclonal antibody targeting human Smad1 (Q15797). Suitable for WB, IHC-P, IP, ELISA; cross-reacts with human and mouse. Ideal for BMP signaling research, developmental biology, and tissue homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Smad1 Rabbit mAb
Remarks/Alias	BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1; d1
Species	Human
Gene ID (Human)	4086
Gene ID	4086
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-240 of human Smad1 (NP_005891.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15797
Protein Weight	52kDa
Shipping	Ice bag