


Anti-Smad1/Smad5/Smad9 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15797/Q99717/O15198

Item Number: CM0024699



Introduction

Rabbit polyclonal antibody against Smad1/Smad5/Smad9, validated for WB, IF/ICC, ELISA in human and mouse samples. Suitable for BMP signaling research, embryonic development, and cell differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	Smad1/Smad5/Smad9 Rabbit pAb
Remarks/Alias	Smad1/Smad5/Smad9
Species	Human
Gene ID (Human)	4086/4090/4093
Gene ID	4086/4090/4093
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 19-125 of human Smad1/Smad5/Smad9(NP_005891.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15797/Q99717/O15198
Protein Weight	52kDa
Shipping	Ice bag