

Anti-Smad2 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, If/Icc, Elisa - Q15796

Item Number: CM0025803



Introduction

KO validated rabbit polyclonal antibody against human Smad2, suitable for WB, IHC-P, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for TGF-beta signaling and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] Smad2 Rabbit pAb
Remarks/Alias	JV18; LDS6; CHTD8; MADH2; MADR2; JV18-1; hMAD-2; hSMAD2; d2
Species	Human
Gene ID (Human)	4087
Gene ID	4087
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Smad2 (NP_005892.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15796
Protein Weight	52kDa
Shipping	Ice bag