


Anti-Smad2/Smad3 Rabbit Monoclonal Antibody For Wb, Ip, Chip, Elisa - Q15796 / P84022

Item Number: CM0011255



Introduction

KO validated anti-SMAD2/SMAD3 rabbit monoclonal antibody for Western blot, immunoprecipitation, chromatin immunoprecipitation, and ELISA. Reacts with human, mouse, rat. Targets TGF-beta/activin signaling transducer SMAD2/SMAD3.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SMAD2/SMAD3 Rabbit mAb
Remarks/Alias	CHTD8; JV18; JV18-1; LDS6; MADH2; MADR2; hMAD-2; hSMAD2; HSPC193; HsT17436; JV15-2; LDS1C; LDS3; MADH3; hMAD-3; hSMAD3; mad3
Species	Human
GeneID (Human)	4087/4088
GeneID	4087/4088
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 198-298 of human SMAD2/SMAD3 (NP_005892.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15796/P84022
Protein Weight	25-52 kDa
Shipping	Ice bag