

Anti-Smad2/Smad3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15796/P84022

Item Number: CM0019884



Introduction

Rabbit polyclonal antibody against SMAD2/SMAD3, suitable for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, rat SMAD2/SMAD3. High specificity for TGF-beta signaling research and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	SMAD2/SMAD3 Rabbit pAb
Remarks/Alias	HSPC193; HsT17436; JV15-2; LDS1C; LDS3; MADH3; SMAD2/SMAD3
Species	Human
GeneID (Human)	4087/4088
GeneID	4087/4088
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 150-270 of human SMAD2/SMAD3 (NP_005893.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15796/P84022
Protein Weight	48kDa/52kDa/25kDa/35kDa/43kDa
Shipping	Ice bag