


Anti-Smad3 (Ko Validated) Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - P84022

Item Number: CM0006476



Introduction

High specificity rabbit monoclonal antibody targeting human Smad3 (P84022). KO validated for Western blot, immunoprecipitation, immunohistochemistry, and ELISA. Suitable for investigating TGF- β -mediated transcription and chondrogenesis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Smad3 Rabbit mAb
Remarks/Alias	LDS3; mad3; LDS1C; MADH3; JV15-2; hMAD-3; hSMAD3; HSPC193; HsT17436; d3
Species	Human
Gene ID (Human)	4088
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 119-219 of human Smad3 (NP_005893.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P84022
Protein Weight	48kDa
Shipping	Ice bag