


Anti-Smad3 Rabbit Monoclonal Antibody (Ko Validated) For Wb, Ihc-P, Elisa - P84022

Item Number: CM0007349



Introduction

Rabbit monoclonal antibody targeting human Smad3, validated for Western blot, immunohistochemistry (paraffin), and ELISA. Detects Smad3 in human, mouse, and rat samples. Ideal for studying TGF- β signaling pathways.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] Smad3 Rabbit mAb
Remarks/Alias	LDS3; mad3; LDS1C; MADH3; JV15-2; hMAD-3; hSMAD3; HSPC193; HsT17436; d3
Species	Human
GeneID (Human)	4088
GeneID	4088
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Smad3 (NP_005893.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84022
Protein Weight	48kDa
Shipping	Ice bag