


Anti-Smad3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P84022

Item Number: CM0022715



Introduction

High-quality rabbit polyclonal anti-Smad3 antibody validated for Western blot, immunofluorescence, and ELISA in human, mouse, and rat. Ideal for TGF-beta signaling and transcriptional regulation studies.

[Learn More](#)

Parameter	Specification
Product Name	Smad3 Rabbit pAb
Remarks/Alias	LDS3, mad3, LDS1C, MADH3, JV15-2, hMAD-3, hSMAD3, HSPC193, HsT17436, Smad3
Species	Human
GeneID (Human)	4088
GeneID	4088
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-250 of human Smad3 (NP_005893.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84022
Protein Weight	48kDa
Shipping	Ice bag