


Anti-Smad5 Rabbit Monoclonal Antibody (Ko Validated) For Wb, Elisa - Q99717

Item Number: CM0007071



Introduction

Anti-Smad5 rabbit monoclonal antibody validated for WB and ELISA, recognizing human, mouse, and rat Smad5. A key transcription factor in BMP signaling, involved in development and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Smad5 Rabbit mAb
Remarks/Alias	DWFC; JV5-1; MADH5; d5
Species	Human
GeneID (Human)	4090
GeneID	4090
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 186-260 of human Smad5 (NP_005894.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99717
Protein Weight	52 kDa
Shipping	Ice bag