


Anti-Smad5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99717

Item Number: CM0023165



Introduction

Rabbit polyclonal antibody against Smad5, validated for WB, IHC-P, ELISA in human, mouse, rat. Recognizes the 52 kDa SMAD family member 5, a key transcriptional regulator in BMP signaling and metabolism.

[Learn More](#)

Parameter	Value
Product Name	Smad5 Rabbit pAb
Remarks/Alias	DWFC; JV5-1; MADH5; Smad5
Species	Human
Gene ID (Human)	4090
Gene ID	4090
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 186-465 of human Smad5 (NP_005894.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99717
Protein Weight	52kDa
Shipping	Ice bag