


Anti-Smad9 Polyclonal Antibody For Ihc-P And Elisa - O15198

Item Number: CM0029579



Introduction

Rabbit polyclonal antibody targeting human SMAD9 (Smad8), verified for IHC-P and ELISA. Detects SMAD9 in human, mouse, and rat samples. Suitable for BMP signaling pathway studies and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	Smad9 Rabbit pAb
Remarks/Alias	PPH2; MADH6; MADH9; SMAD8; SMAD8A; SMAD8B; SMAD8/9; Smad9
Species	Human
GeneID (Human)	4093
GeneID	4093
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Smad9 (NP_005896.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15198
Protein Weight	52kDa
Shipping	Ice bag