


Anti-Smap1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Iyb5

Item Number: CM0021449



Introduction

High-quality anti-SMAP1 rabbit polyclonal antibody for Western blot, immunofluorescence, and ELISA. Cross-reacts with human, mouse, and rat SMAP1. Ideal for studying clathrin-dependent endocytosis and erythropoiesis.

[Learn More](#)

Parameter	Value
Product Name	SMAP1 Rabbit pAb
Remarks/Alias	SMAP-1; SMAP1
Species	Human
Gene ID	60682
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IYB5
Protein Weight	50kDa
Shipping	Ice bag