

Anti-Smarca1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P28370

Item Number: CM0013179

Introduction

Rabbit polyclonal antibody against human SMARCA1 (SWI/SNF-related chromatin remodeler). Validated for WB, IHC-P, IF/ICC, IP, ELISA. Detects human, mouse, rat SMARCA1. Ideal for epigenetic and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	SMARCA1 Rabbit pAb
Remarks/Alias	SWI; ISWI; SWI2; SNF2L; SNF2L1; SNF2LB; SNF2LT; hSNF2L; NURF140; SMARCA1
Species	Human
GeneID (human)	6594
GeneID	6594
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human SMARCA1 (NP_003060.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28370
Protein Weight	123kDa
Shipping	Ice bag