

Anti-Smarca2 / Brm Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Chip, Elisa - P51531

Item Number: CM0007054



Introduction

Rabbit monoclonal antibody against human SMARCA2/BRM, validated for WB, IHC-P, IF/ICC, IP, ChIP, ELISA. Target ATPase involved in chromatin remodeling, useful for transcription and neurodevelopment research.

[Learn More](#)

Parameter	Value
Product Name	SMARCA2 / BRM Rabbit mAb
Remarks/Alias	BIS; BRM; SNF2; SWI2; hBRM; NCBRS; Sth1p; BAF190; SNF2L2; SNF2LA; hSNF2a; SMARCA2 / BRM
Species	Human
GeneID (Human)	6595
GeneID	6595
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ChIP, ELISA
Cross Reactivity	Human
SWISS	P51531
Protein Weight	181kDa
Shipping	Ice bag