

Anti-Smarca5/Snf2H Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O60264

Item Number: CM0010255



Introduction

High-quality rabbit monoclonal antibody against human SMARCA5/SNF2H (O60264). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for chromatin remodeling and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	SMARCA5/SNF2H Rabbit mAb
Remarks/Alias	ISWI, SNF2H, hISWI, hSNF2H, WCRF135, SMARCA5/SNF2H
Species	Human
GeneID (Human)	8467
GeneID	8467
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 953-1052 of human SMARCA5/SNF2H (O60264).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60264
Protein Weight	122kDa
Shipping	Ice bag