

Anti-Smarca5/Snf2H Rabbit Polyclonal Antibody For Wb And Elisa - O60264

Item Number: CM0020326



Introduction

Rabbit polyclonal antibody targeting human SMARCA5/SNF2H, an ATP-dependent chromatin remodeler. Validated for WB and ELISA; cross-reacts with mouse. Ideal for studying nucleosome remodeling, DNA replication, and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	SMARCA5/SNF2H Rabbit pAb
Remarks/Alias	ISWI; SNF2H; hISWI; hSNF2H; WCRF135; SMARCA5/SNF2H
Species	Human
Gene ID (Human)	8467
Gene ID	8467
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-400 of human SMARCA5/SNF2H (NP_003592.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60264
Protein Weight	122kDa
Shipping	Ice bag