


Anti-Smarcad1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H4L7

Item Number: CM0028774



Introduction

SMARCAD1 Rabbit pAb validated for WB, IF/ICC, ELISA. Targets human, mouse, rat. Ideal for chromatin remodeling and DNA repair research. 117 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	SMARCAD1 Rabbit pAb
Remarks/Alias	HRZ; TYS; ETL1; HEL1; ADERM; BASNS; SMARCAD1
Species	Human
Gene ID (Human)	56916
Gene ID	56916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SMARCAD1 (NP_001121901.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4L7
Protein Weight	117kDa
Shipping	Ice bag