

Anti-Smarcb1/Snf5 Polyclonal Antibody For Wb, Elisa - Q12824

Item Number: CM0022078



Introduction

SMARCB1/SNF5 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with Human, Mouse, and Rat. Core component of the BAF chromatin remodeling complex with roles in cell cycle control and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	SMARCB1/SNF5 Rabbit pAb
Remarks/Alias	RDT; CSS3; INI1; SNF5; Snr1; BAF47; INI-1; MRD15; RTPS1; Sfh1p; hSNF5; SNF5L1; SWNTS1; PPP1R144; SMARCB1/SNF5
Species	Human
GeneID (Human)	6598
GeneID	6598
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12824
Protein Weight	44kDa
Shipping	Ice bag