

Anti-Smarcc1/Baf155 Monoclonal Antibody For Wb, If/Icc, Elisa, Chip - Q92922

Item Number: CM0011419



Introduction

KO Validated recombinant rabbit monoclonal antibody targeting SMARCC1/BAF155, suitable for WB, IF/ICC, ELISA, and ChIP. Detects human and mouse SMARCC1 (~123 kDa). Essential for chromatin remodeling and neural development.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SMARCC1/BAF155 Rabbit mAb
Remarks/Alias	HYC5; Rsc8; SRG3; SWI3; BAF155; CRACC1; 55
Species	Human
GeneID (human)	6599
GeneID	6599
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse
SWISS	Q92922
Protein Weight	123kDa
Shipping	Ice bag