

Anti-Smarcc2 Polyclonal Antibody For Wb, Elisa - Q8Taq2

Item Number: CM0022498



Introduction

Rabbit polyclonal antibody against SMARCC2 (BAF170), a SWI/SNF chromatin remodeling subunit. Suitable for WB and ELISA. Detects human, mouse, and rat SMARCC2. Supports research into transcriptional regulation, neurogenesis, and myeloid differentiation.

[Learn More](#)

Parameter	Value
Product Name	SMARCC2 Rabbit pAb
Remarks/Alias	CSS8, Rsc8, BAF170, CRACC2, SMARCC2
Species	Human
GeneID	6601
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human SMARCC2 (NP_620706.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TAQ2
Protein Weight	133kDa
Shipping	Ice bag