

Anti-Smarcd1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Gm5

Item Number: CM0024887



Introduction

Rabbit polyclonal antibody targeting SMARCD1 (BAF60A), a SWI/SNF chromatin remodeling subunit; confirmed for WB and ELISA in human, mouse, and rat samples. Ideal for transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	SMARCD1 Rabbit pAb
Remarks/Alias	CSS11; Rsc6p; BAF60A; CRACD1; SMARCD1
Species	Human
GeneID (Human)	6602
GeneID	6602
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 216-515 of human SMARCD1 (NP_003067.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96GM5
Protein Weight	58kDa
Shipping	Ice bag