

Anti-Smarcd2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92925

Item Number: CM0015805

Introduction

Rabbit polyclonal antibody targeting SMARCD2 (BAF60B), a chromatin remodeling factor. Validated for WB, IHC-P, and ELISA in human, mouse, and rat. Ideal for chromatin and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	SMARCD2 Rabbit pAb
Remarks/Alias	SGD2; Rsc6p; BAF60B; CRACD2; PRO2451; SMARCD2
Species	Human
GeneID (Human)	6603
GeneID	6603
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 454-531 of human SMARCD2 (NP_001091896.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92925
Protein Weight	59kDa
Shipping	Ice bag