

Anti-Smarcd3 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Ste5

Item Number: CM0024518

Introduction

High-quality rabbit polyclonal antibody targeting SMARCD3 (BAF60C), a key chromatin remodeling factor. Validated for WB and ELISA, reacts with human, mouse, and rat. Ideal for transcription regulation and neurogenesis research.

[Learn More](#)

Parameter	Value
Product Name	SMARCD3 Rabbit pAb
Remarks/Alias	Rsc6p; BAF60C; CRACD3; SMARCD3
Species	Human
GeneID (Human)	6604
GeneID	6604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 275-440 of human SMARCD3 (NP_001003801.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6STE5
Protein Weight	55kDa
Shipping	Ice bag