


Anti-Smc1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q14683

Item Number: CM0008870



Introduction

High-quality rabbit monoclonal anti-SMC1 antibody validated for Western blot, IF/ICC, and ELISA. Reacts with human and mouse SMC1 (Q14683). Ideal for chromosome cohesion and DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	SMC1 Rabbit mAb
Remarks/Alias	SMC1; SMCB; CDLS2; DEE85; SB1.8; EIEE85; SMC1L1; DXS423E; SMC1alpha
Species	Human
GeneID (Human)	8243
GeneID	8243
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1134-1233 of human SMC1 (Q14683).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14683
Protein Weight	143kDa
Shipping	Ice bag