


## Anti-Smc1 Rabbit Polyclonal Antibody For Wb And Elisa - Q14683

Item Number: CM0025015



### Introduction

Rabbit polyclonal antibody against human SMC1, validated for WB and ELISA. Detects endogenous SMC1 protein (~143 kDa) across human, mouse, and rat samples. Suitable for chromosome cohesion and DNA repair research.

[Learn More](#)

| Parameter        | Value                                                                                                                                |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | SMC1 Rabbit pAb                                                                                                                      |
| Remarks/Alias    | SMC1; SMCB; CDLS2; DEE85; SB1.8; EIEE85; SMC1L1; DXS423E; SMC1alpha                                                                  |
| Species          | Human                                                                                                                                |
| GeneID (Human)   | 8243                                                                                                                                 |
| GeneID           | 8243                                                                                                                                 |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human SMC1 (NP_006297.2). |
| Source           | Rabbit                                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                                |
| Application      | WB, ELISA                                                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                    |
| SWISS            | Q14683                                                                                                                               |
| Protein Weight   | 143kDa                                                                                                                               |
| Shipping         | Ice bag                                                                                                                              |