


Anti-Smc3 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q9Uqe7

Item Number: CM0002059



Introduction

High-quality rabbit monoclonal SMC3 antibody validated for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Suitable for cohesin complex and chromosome cohesion studies.

[Learn More](#)

Parameter	Value
Product Name	SMC3 Rabbit mAb
Remarks/Alias	BAM; BMH; HCAP; CDLS3; CSPG6; SMC3L1; SMC3
Species	Human
Gene ID (Human)	9126
Gene ID	9126
Immunogen	Synthetic Peptide / A synthetic peptide corresponding to a sequence within amino acids 1118-1217 of human SMC3 (NP_005436.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQE7
Protein Weight	142kDa
Shipping	Ice bag