


## Anti-Smc3 Rabbit Polyclonal Antibody (Smc3 Pab) For Wb, Ip, Elisa - Q9Uqe7

Item Number: CM0023246



### Introduction

Rabbit polyclonal antibody targeting human SMC3, validated for WB, IP, and ELISA. Reacts with human and mouse SMC3. Suitable for chromosome cohesion, cell cycle, and DNA repair research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | SMC3 Rabbit pAb                                                                                                    |
| Remarks/Alias    | BAM, BMH, HCAP, CDLS3, CSPG6, SMC3L1, SMC3                                                                         |
| Species          | Human                                                                                                              |
| GeneID           | 9126                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 909-1199 of human SMC3 (NP_005436.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IP, ELISA                                                                                                      |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | Q9UQE7                                                                                                             |
| Protein Weight   | 142kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                            |