


Anti-Smc4 Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Ntj3

Item Number: CM0004714



Introduction

Anti-SMC4 Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA; reacts with Human, Mouse, Rat; detects condensin complex protein SMC4, ~147kDa. Ideal for chromosome biology research.

[Learn More](#)

Parameter	Value
Product Name	SMC4 Rabbit mAb
Remarks/Alias	CAPC; CAP-C; SMC-4; SMC4L1; SMC4
Species	Human
GeneID (Human)	10051
GeneID	10051
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1189-1288 of human SMC4 (Q9NTJ3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NTJ3
Protein Weight	147kDa
Shipping	Ice bag