

Anti-Smc5 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Iy18

Item Number: CM0019391

Introduction

Rabbit polyclonal antibody targeting human SMC5, validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for DNA repair and telomere biology research.

[Learn More](#)

Parameter	Value
Product Name	SMC5 Rabbit pAb
Remarks/Alias	ATELS2; SMC5L1; SMC5
Species	Human
Gene ID (Human)	23137
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 872-1101 of human SMC5 (NP_055925.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IY18
Protein Weight	129kDa
Shipping	Ice bag