


Anti-Smchd1 Polyclonal Antibody For Wb, If/Icc, Elisa - A6Nhr9

Item Number: CM0025626



Introduction

Rabbit polyclonal antibody against SMCHD1, suitable for WB, IF/ICC, and ELISA. Recognizes human and rat SMCHD1 (226 kDa). Ideal for chromatin silencing, X-inactivation, and DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	SMCHD1 Rabbit pAb
Remarks/Alias	BAMS; FSHD2; SMCHD1
Species	Human
Gene ID (Human)	23347
Gene ID	23347
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1756-2005 of human SMCHD1 (NP_056110.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	A6NHR9
Protein Weight	226kDa
Shipping	Ice bag