


Anti-Smg1 Rabbit Monoclonal Antibody For Wb And Elisa - Q96Q15

Item Number: CM0002672



Introduction

Rabbit monoclonal antibody targeting human SMG1, validated for Western blot and ELISA. Cross-reacts with mouse. Ideal for studying nonsense-mediated decay and genotoxic stress pathways.

[Learn More](#)

Parameter	Value
Product Name	SMG1 Rabbit mAb
Remarks/Alias	ATX; LIP; 61E3.4; SMG1
Species	Human
GeneID (Human)	23049
GeneID	23049
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 3562-3661 of human SMG1 (Q96Q15).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96Q15
Protein Weight	411 kDa
Shipping	Ice bag