


Anti-Smg7 Rabbit Polyclonal Antibody For Wb, Elisa - Q92540

Item Number: CM0018666



Introduction

Rabbit polyclonal antibody against SMG7, a key nonsense-mediated mRNA decay factor. Validated for WB and ELISA, it cross-reacts with human, mouse, and rat samples. Ideal for mRNA surveillance and post-transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	SMG7 Rabbit pAb
Remarks/Alias	EST1C, SGA56M, C1orf16, SMG7
Species	Human
GeneID	9887
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SMG7 (NP_963863.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92540
Protein Weight	127kDa
Shipping	Ice bag