


Anti-Smmhc/Myh11 Rabbit Polyclonal Antibody For Wb - P35749

Item Number: CM0016806



Introduction

Rabbit polyclonal antibody targeting human SMMHC/MYH11, validated in WB, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for smooth muscle contraction research.

[Learn More](#)

Parameter	Value
Product Name	SMMHC/MYH11 Rabbit pAb
Remarks/Alias	AAT4; FAA4; SMHC; SMMHC; VSCM2; SMMS-1; SMMHC/MYH11
Species	Human
GeneID (Human)	4629
GeneID	4629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1680-1979 of human SMMHC/SMMHC/MYH11 (NP_001035203.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35749
Protein Weight	227kDa
Shipping	Ice bag