


Anti-Smn1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q16637

Item Number: CM0010728



Introduction

Rabbit monoclonal antibody against human SMN1, validated for Western blot, IHC-P, and ELISA. Cross-reacts with rat. Suitable for studying spliceosome assembly, snRNP biogenesis, and SMA-related research.

[Learn More](#)

Parameter	Value
Product Name	SMN1 Rabbit mAb
Remarks/Alias	SMA; SMN; SMA1; SMA2; SMA3; SMA4; SMA@; SMNT; BCD541; GEMIN1; TDRD16A; T-BCD541; SMN1
Species	Human
Gene ID (Human)	6606
Gene ID	6606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 103-200 of human SMN1 (Q16637).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q16637
Protein Weight	32kDa
Shipping	Ice bag