


Anti-Smn2 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q16637

Item Number: CM0015353



Introduction

Rabbit polyclonal antibody against human SMN2, suitable for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for studying spliceosomal snRNP assembly and motor neuron development.

[Learn More](#)

Parameter	Value
Product Name	SMN2 Rabbit pAb
Remarks/Alias	SMNC; BCD541; GEMIN1; TDRD16B; C-BCD541; SMN2
Species	Human
Gene ID	6607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human SMN2 (NP_059107.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16637
Protein Weight	32kDa
Shipping	Ice bag