

Anti-Smn2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16637

Item Number: CM0012841



Introduction

Rabbit polyclonal antibody targeting human SMN2 (survival motor neuron protein). Suitable for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for investigating splicing, neuromuscular disorders, and spinal muscular atrophy mechanisms.

[Learn More](#)

Parameter	Value
Product Name	SMN2 Rabbit pAb
Remarks/Alias	SMNC, BCD541, GEMIN1, TDRD16B, C-BCD541, SMN2
Species	Human
GeneID (Human)	6607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human SMN2 (NP_059107.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16637
Protein Weight	32kDa
Shipping	Ice bag