

Anti-Smoc1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H4F8

Item Number: CM0023163

Introduction

SMOC1 Rabbit pAb is a polyclonal antibody that recognizes human and mouse SMOC1 protein. Suitable for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. SMOC1 is a secreted calcium-binding protein essential for eye and limb development.

[Learn More](#)

Parameter	Value
Product Name	SMOC1 Rabbit pAb
Remarks/Alias	OAS; SMOC1
Species	Human
GeneID (Human)	64093
GeneID	64093
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-250 of human SMOC1 (NP_071420.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H4F8
Protein Weight	48kDa
Shipping	Ice bag